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## Customer Reprocessing Notice Regarding Mildly Alkaline Cleaning Agents

December 4, 2012

SUBJECT: Use of Mildly Alkaline Cleaning Agents (pH  $\leq$  11) During the Reprocessing of Reusable Instruments used with *da Vinci*<sup>®</sup> Surgical Systems

Neu: „Mildly Alkaline Enzymatic Cleaning Agents“

Dear Customer,

In response to requests from customers, we have performed testing to determine if mildly alkaline cleaning agents (defined as alkaline agents with a pH  $\leq$  11) can be used in the reprocessing of 8mm reusable instruments used with the *da Vinci* Surgical system. Our reprocessing instructions currently recommend the use of pH-neutral cleaning agents only.

Our testing showed that 8mm reusable instruments are compatible with mildly alkaline cleaners. Our tests showed no decrease in instrument performance or number of lives when mildly alkaline cleaners are used. Alkaline cleaners without silicates may cause cosmetic surface roughening and discoloration of the aluminum parts found inside the back-end of the instrument; however, this did not affect instrument performance or the number of lives.

Our testing relates to reusable 8mm instruments only. It does not include 5mm instruments, the EndoWrist<sup>®</sup> Stabilizer or reusable accessories for the Standard, S or Si *da Vinci* Systems.

We are in the process of assessing changes to our reprocessing instructions (Instructions for Use) as a result of these tests, but we wanted to respond to customers with these test results.

If you have questions regarding this message, please contact Customer Service at the appropriate number below.

US: 1-800-876-1310

EU: +800.0821.2020 or +41.21.821.2020